

A New Method for the Quantitation of a Large Number of Pesticides Residues in Food by Gas Chromatography Mass Spectrometry

V. C. Paul¹, S. Prithamson¹, S. L. Sengar¹, M. D. Mishra², D. Wang²

¹ Statutory Analytical Instrumentation, IARI, New Delhi, India

² Scientific Analytical Instrumentation, IARI, New Delhi, India

Scientific Analytical Instrumentation, IARI, New Delhi, India

Scientific Analytical Instrumentation, IARI, New Delhi, India

Introduction

There is a requirement within the European Union (EU) and elsewhere to measure the levels of certain pesticide residues in food. The United Kingdom government has issued a series of statutory instruments that specify maximum pesticide residue levels in crops, food and animal feeding stuffs. This legislation incorporates FAO/WHO CODEX Alimentarius commission methods of analysis and sampling. In practice some 62 pesticide residues are routinely measured down to 100 μg of lot weight, which is normally carried out by a number of separate specific GC-MS methods. The statistical analysis of the results from these methods using the GC-MS method that determines all the target compounds and allows semi-quantitative analysis of other residues that may be present.

Method

The analysis was performed on extracts taken from homogenised feedstuffs. The extraction method was based on the FDA Pesticide Analytical Manual section 302 and consists of taking a 100g representative sample which was spiked with an internal standard mixture at a level of 0.1 mg/kg and then extracted at room temperature with 100ml of acetone. This mixture was then filtered and the acetone/water phase back extracted with 200ml of 50:50 Dichloroethane ether. The liquid/liquid extraction was then repeated a further two times with DCM. The resulting combined extracts were then analysed on a GC-MS system. The GC is connected to a 30m x 0.25mm DB-5ms. The GC is connected to a VG 70S mass spectrometer fitted with an electron ionisation source. Full mass spectra were collected and recorded over a mass range of 65–450 Daltons.

Results

The results show that a wide range of pesticide residues can be analysed in a single GC-MS run, avoiding the need to carry out multiple injections on instruments with selective detectors. Table 1 shows the list of targeted pesticides and internal standards used for the analysis. In Figure 1 we can see the complexity of the standard mixture. Figure 2 gives a good illustration of the power of this technique on a grape extract. Note the peak at a retention time of 26.5 min, identified as pyridoxine, which is present well below the limit of 5 mg/kg and quantified at 0.02 mg/kg. Table 2 shows the dynamic range of the analysis is clear. The spectrum for the target is shown in Figure 3 and indicates the dynamic range of the system given the excellent spectral quality of this component at this level. A second sample of grapes shows a component Figure 4 that is not in the target database but is a pesticide. curball, the component's spectrum and its library match are shown in Figure 5.

No	Compound	Peak (min)	MW	Library	No	Compound	Peak (min)	MW	Library
1	1,1-Dichloroethane	1.8	98	1,1-Dichloroethane	31	Phenylpyrazole	26.5	149	Phenylpyrazole
2	Acetone	2.1	58	Acetone	32	Pyridoxine	26.5	153	Pyridoxine
3	Acetone	2.1	58	Acetone	33	Acetone	2.1	58	Acetone
4	Acetone	2.1	58	Acetone	34	Acetone	2.1	58	Acetone
5	Acetone	2.1	58	Acetone	35	Acetone	2.1	58	Acetone
6	Acetone	2.1	58	Acetone	36	Acetone	2.1	58	Acetone
7	Acetone	2.1	58	Acetone	37	Acetone	2.1	58	Acetone
8	Acetone	2.1	58	Acetone	38	Acetone	2.1	58	Acetone
9	Acetone	2.1	58	Acetone	39	Acetone	2.1	58	Acetone
10	Acetone	2.1	58	Acetone	40	Acetone	2.1	58	Acetone
11	Acetone	2.1	58	Acetone	41	Acetone	2.1	58	Acetone
12	Acetone	2.1	58	Acetone	42	Acetone	2.1	58	Acetone
13	Acetone	2.1	58	Acetone	43	Acetone	2.1	58	Acetone
14	Acetone	2.1	58	Acetone	44	Acetone	2.1	58	Acetone
15	Acetone	2.1	58	Acetone	45	Acetone	2.1	58	Acetone
16	Acetone	2.1	58	Acetone	46	Acetone	2.1	58	Acetone
17	Acetone	2.1	58	Acetone	47	Acetone	2.1	58	Acetone
18	Acetone	2.1	58	Acetone	48	Acetone	2.1	58	Acetone
19	Acetone	2.1	58	Acetone	49	Acetone	2.1	58	Acetone
20	Acetone	2.1	58	Acetone	50	Acetone	2.1	58	Acetone
21	Acetone	2.1	58	Acetone	51	Acetone	2.1	58	Acetone
22	Acetone	2.1	58	Acetone	52	Acetone	2.1	58	Acetone
23	Acetone	2.1	58	Acetone	53	Acetone	2.1	58	Acetone
24	Acetone	2.1	58	Acetone	54	Acetone	2.1	58	Acetone
25	Acetone	2.1	58	Acetone	55	Acetone	2.1	58	Acetone
26	Acetone	2.1	58	Acetone	56	Acetone	2.1	58	Acetone
27	Acetone	2.1	58	Acetone	57	Acetone	2.1	58	Acetone
28	Acetone	2.1	58	Acetone	58	Acetone	2.1	58	Acetone
29	Acetone	2.1	58	Acetone	59	Acetone	2.1	58	Acetone
30	Acetone	2.1	58	Acetone	60	Acetone	2.1	58	Acetone

Table 1 List of target compounds and internal standards




Figure 1: Chromatogram of a standard mixture




Figure 2: Chromatogram of a grape sample showing low level




Figure 3: Spectrum of a component from a grape




Figure 4: Chromatogram of a grape sample showing low level




Figure 5: Mass spectrum of curball, a component not in the target database

Conclusion

This method of screening for target pesticides has been shown to both rapid and sensitive when compared to other published methods using a range of selective detectors. The use of a magnetic sector mass spectrometer for this application over a lower cost quadrupole mass spectrometer provides a number of advantages. The ability to collect full spectra is essential if targets not in the list are to be reported. The authors feel that this application would be well suited to high-throughput instrumentation that would be capable of reducing full mass spectra whilst not compromising sensitivity. There is great interest in reducing GC run times using new multiple GC columns, an advance that cannot be utilised using TOF instruments. The dynamic range of such a system must however be at least equivalent to the current instrumentation.

References

1. CODEX Alimentarius Volume 2A Part 1
2. Statutory Instrument 1989 NO 2483 Pesticide Residue Levels in Crops, food and feeding stuffs