

A Comparison of N-terminal Sequencing Using Different Sulphonation Reagents with MALDI-PSD-TOFMS

David J Evason, Mark D Mills, Victor C Parr, Alexis J Polley : Scientific Analysis Instruments, Manchester, UK.

Overview

Charge directed fragmentation via N-terminal sulphonation

- Fragmentation occurs more readily
- Predictable ladder sequence (y ion only, for tryptics)
- Via a "harmonic field reflectron" MALDI-TOFMS

APPROACHES

- 4-sulphophenyl isothiocyanate (SPITC)
- 2-sulphophenyl benzoic acid cyclic anhydride (SBAn)
- 2-sulfopropionic acid NHS-ester (CAF)

CONCLUSIONS

- SBAn is difficult to use, derivatives aren't seen in reflectron
- SPITC had the best performance, is more robust than CAF and is cheaper

Introduction

Charge directed fragmentation (CDF) increases the fragmentation rate of peptides and directs their fragmentation into a predictable pattern. For this purpose, sulphonation has proved useful. The following presentation compares three methods of sulphonation for the purpose of peptide sequencing by a harmonic field reflectron MALDI-PSD-TOFMS.

Methods

All materials were from Sigma, Poole, UK, except for CAF reagent (gift from M Liming, Amersham Biosciences).

Protein digested with 5% w/v DPC treated trypsin, in 0.05M ammonium bicarbonate (AMBC).

Sulphonation a) *SPITC* : sample diluted 1:1 with 50mg/ml 4-sulphophenyl isothiocyanate in 0.04M AMBC. b) *SBAn*: sample diluted 1:1 with 10mg/ml 2-sulfolobenzoic acid cyclic anhydride in dry acetone. c) *CAF*: sample diluted 1:1 in 10mg/ml CAF reagent in 0.25M Na₂CO₃ pH 9.4. Samples incubated for 2 hours at room temperature. Sample clean up prior to analysis on C-18 z-tips (Millipore).

All MALDI-TOFMS was performed on a *LaserToF RS*, SAI, UK, using a high performance, harmonic field reflectron MALDI-TOFMS (Thompson et al).

Results




Figure 1: Yield of derivatization in linear




Figure 2: Yield of derivatization in reflectron




Figure 3: PSD of (top left, clockwise) 160fmol, FibA, 80fmol SPITC-FibA, 16fmol FibA, 8fmol SPITC-FibA




Figure 4: PSD of (top left, clockwise) 80fmol SBAn-FibA, 80fmol SPA-FibA, 8fmol SBAn-FibA, 8fmol SPA-FibA

	1	2	3	4	5
None	77.9	14, 11, 9, 7 and 1	na	na	Sig
SBAn	78.3	15, 14, 13 and 11	40	0.02	Sig
SPITC	100.0	All	53	8.1	Y 15
CAF	100.0	15 only	43	0.49	none

Table 1: Comparison of Sulphonation Treatments on PSD Efficacy

- Percentage of ions that are C-terminal in PSD of 8 pmol peptide
- Y ions present at the lowest dilution tested. (8 fmol for all derivatizations, 16 fmol for underivatized)
- Ratio of derivatized to underivatized linear mode
- Ratio of derivatized to underivatized reflectron mode
- ISD - In Source Decay. Sig = significant

Discussion

Figures 1 and 2 compare the sulphonation chemistries in linear and reflectron modes. The linear mode cannot distinguish between a derivative ion and its in-flight fragments (metastable ions). It therefore gives a better picture of yield since the derivative peak will be composed of the derivative ions at the time of their ionisation in the source. The reflectron spectra is more complicated. It will have peaks of the derivative, and should it fragment during its flight before the reflectron, it will have peaks of these fragments (metastable ions). The reflectron spectra also shows that sulphonation has caused the molecules to fragment in the source (ISD ions). Timing is everything. To be of maximum use in PSD, the derivative must fragment, during its flight to the reflectron, completely, into an equal distribution of all possible y-ions. From the reflectron spectra it appears that CAF and SBAn have done this. However, both have poorer distributions and yields of y-ions at low peptide concentrations (Figure 4). For SBAn derivatives, this is because fragmentation occurs too early, in the source (ISD ions Fig 2).

Conclusion

A good sulphonation (CDF) agent should: give y-only spectra, give readily-identifiable derivatives (for automation), have good y-ion sequences at low concentrations of peptide. In this study SPITC was found to be the best sulphonation agent for PSD. It gave a good yield of readily identifiable derivatives (Fig 1 & 2), had little in source decay (Fig 2), gave a continuous and full y-ion sequence at 80 fmol (Fig 3) and it gave a continuous y-ion sequence from y-15 to y-4 at 8 fmol (Fig 3).

As a demonstration SPITC and SBAn derivatization of a cytochrome c tryptic peptide are shown (Fig 5). The native peptide had less useful PSD (data not shown).

References

- Chen et al Rapid Comms in MS 2004: 18,191-198.
- Gevaert et al 2001 Electrophoresis 22: 1645
- Keough et al 1999 PNAS 96: 7131
- Marekov and Steimert J Mass Spectrom. 2003, 38: 373-377
- Thompson et al, Proc. 47th ASMS, 1999, Dallas Tx, p 3066
- Wang, Kalb and Cotter Rapid Comms in MS 2004: 18,96-102




Figure 5: CDF of a tryptic peptide from digested cytochrome c. SPITC v SBAn